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
Titripur® - The Measure of all Things - Titration Guide for Titer Determination

2.2 Sodium thiosulfate solution 0.1 mol/l (0.1 N) (Catalogue No. 109147)

Reaction type: Redox titration/iodometry
Kind of titration: Back-titration based on substitution
Volumetric standard: Potassium iodate
 Potassium iodate should be preferred to the toxic potassium dichromate which is suspected to be carcinogenic.
Method of initial introduction: Division into aliquots by gravimetry

Indication

a) **visual:** with starch solution 1 %; colour change from blue to colourless,
 b) **potentiometrical:** with platinum single rod measuring chain and automatic titrator.

Analytical procedure

a) **Preparation of a potassium iodate solution 1/60 mol/l (0.1 N) (Catalogue No. 102404)**
 Approx. 3.57 g of potassium iodate (Volumetric Standard), dried for 2 h at 130 °C, are exactly weighed, transferred to a 1-l volumetric flask, dissolved in water GR for analysis and made up to 1 l with water. The total weight is determined to the second decimal.

Calculation of the mass fraction w (KIO₃) in mg/g solution:

$$w \text{ (mg KIO}_3\text{/g solution)} = \frac{\text{weight of volumetric standard (mg)} \cdot Cf}{\text{total weight (g)}}$$

Cf = correction factor, calculated from the effective content of the volumetric standard acc. to the certificate divided by 100.

Volumetric-methode	Standard-Reference-material (SRM)	Volumetric Solution	Volumetric standards Merck
Acidimetry	Tris(hydroxymethyl)-aminomethane (NIST)	Hydrochloric acid	Tris(hydroxymethyl)-aminomethane (NIST) Order No. 1.02408.0080
		Hydrochloric acid	Sodium carbonate Order No. 1.02405.0080
Alkalimetry	Potassium hydrogen-phthalate (NIST)	Sodium hydroxid solution	Potassium hydrogen phthalat Order No. 1.02400.0080
	Benzoic acid (NIST)	Sodium hydroxid solution	Benzoic Acid Order No. 1.02401.0060
Argentometry	Potassium chloride (NIST)	Silver nitrate solution	Sodium chloride Order No. 1.02406.0080
Reductometry	Potassium dichromate (NIST)	Sodium thiosulfate solution	Potassium dichromate Order No. 1.02403.0100
		Sodium thiosulfate solution	Potassium iodate Order No. 1.02404.0100
Oxidimetry	di-Sodium oxalate (NIST)	Potassium permanganate solution	Fe(II)-ethylen-diammonium sulfate Order No. 1.02402.0080
		Potassium permanganate solution	di-Sodium oxalate Order No. 1.02407.0060
Complexometry	Zinc (NIST)	EDTA (Titriplex® III-solution)	Zinc (NIST) Order No. 1.02409.0100
		EDTA (Titriplex® III-solution)	Calcium carbonate Order No. 1.02410.0100

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
Stability of sodium hydroxide volumetric solution (item number: 1.09137) in a Titripac®:

Volumetric solutions contain reagents that react quantitatively with the substance to be measured in the sample solution. The concentration of the volumetric solution is indicated by the molarity (mol reagent per litre of solution). For the optimal consumption per titration the concentration of the volumetric solution should be close according to the expected concentration of the sample. A standardization (titer determination) with a volumetric standard is recommended to correct for all influences on the measurement result arising from lab conditions, equipment and handling.

To demonstrate the stability of this sodium hydroxide volumetric solution the Titripac® was tested under daily routine conditions.

Every six months the concentration of the sodium hydroxide solution was determined by titration with hydrochloric acid standard solution (item number 1.09057) which was standardized against volumetric standard triiodoacetic acid (item number 1.02408).

The following diagrams show the measured concentration of sodium hydroxide solution (item number 1.09137) in a Titripac 4L and 10L over the shelf life of 36 months.

Ordering information:

Sodium hydroxide solution
Sodium hydroxide solution

Ordering Number:

1.09137.4000 (4L-Titripac®)
1.09137.9010 (10L-Titripac®)

Technical Data Sheet

The diagrams on page 1 show the concentration of sodium hydroxide volumetric solution 1 mol/l in between the range of 0.999 mol/l and 1.001 mol/l over the shelf life of three years. The solution is stable in both packages (Titripac 4L and 10L). The calculation of the concentration is done according to the following equation:

$$c(x)_{\text{actual value}} = c(x)_{\text{nominal value}} \cdot z$$

$c(x)$ indicated value = concentration as specified by the supplier
 z = titer of the volumetric solution

$c(x)$ actual value = concentration after correction for the titer

The titer as correction factor depends on lab conditions, equipment and handling.

The titer determination is carried out with volumetric standards or with standard solutions calibrated directly against volumetric standards.

Volumetric solutions are temperature sensitive. Therefore the titer changes with the temperature. Thus it is recommended to do the titer determination at the same temperature as the following sample titration.

The titer determination of a volumetric solution is essential for precise sample measurements. Errors due to titer determination affect the following measurements as systematic deviation.

Instructions for use of a Titripac®:

Open the Titripac® by pressing on the pre-cutted area.

(Do not use a knife for opening because you can easily damage the liner bag of the Titripac®.) Take the installed tap out of the box and fix it by closing the pre-cutted area.

By opening the tap, solution can be withdrawn without the risk of contamination.

Finally the Titripac® can be connected directly to the titrator by means of a small adapter and a hose with a thread on both sides.

Hose and adapter can be ordered under item number 1.88075.0001.

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Titripur®
Reag. Ph Eur, Reag. USP

Perchloric acid

in anhydrous acetic acid, for titrations in non-aqueous media
(contains acetic acid)

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Ácido perclórico

(contiene Ácido acético)

Acide perchlorique

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Ácido perclórico

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Acido perclorico

(contiene Acido acetico)

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
c(HClO₄) = 0.1 mol/l (0.1 N)

Made in Germany

Merck KGaA, 64271 Darmstadt
Germany, Tel. +49(0)6151 72-2440
EMD Millipore Corporation
400 Summit Drive, Burlington MA 01803
USA, Tel. +1-978-715-4321
sigmaaldrich.com

4 022536 110295

© 2018

Danger. Flammable liquid and vapour. May be corrosive to metals. Causes severe skin burns and eye damage. Keep away from heat. Wear protective gloves/ protective clothing/ eye protection/ face protection. IF SWALLOWED: Rinse mouth. Do NOT induce vomiting. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. IF exposed or concerned: immediately call a POISON CENTER or doctor/ physician. | Adjusted to titer t = 1,000 at 20 °C. Variation: approx. 0.1 %/°C. Cannot be stored indefinitely. The solution crystallizes at below 18°C. The crystals melt reversibly upon heating.

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