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
Application Note MK_AN6557EN Ver. 1.0

Application Note

Impact of the new EN 17141 on the ISO 14698 Validation of MAS-100[®] Air Samplers

General Requirements and Validation

The scope of the new norm is to establish the requirements, recommendations and methodology for microbiological contamination control in clean controlled environments. The document is limited to viable microbial contamination and refers to ISO 14644 about non-viable contamination. ISO 14698 is valid throughout the rest of the world while EN 17141 2020 will replace it in the European Union.

The basis for an efficient microbial air sampler is a good physical and biological efficacy:

"The physical collection efficiency is the cut-off size (d₅₀ value) which defines the aerodynamic equivalent particle diameter size at which the sampler collects 50 % of the particles in the air [12] [13]. The average equivalent diameters of microbe-carrying particles (MCPs) that form the cfu are generally larger than 1

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Table 1. Comparison of requirements for sampling devices and microbial air samplers in EN 17141 and ISO 14698:

Requirements EN 17141	Corresponding Requirements out of ISO 14698
Accessibility into the clean controlled environment for the sampling device	<ul style="list-style-type: none"> ♦ Ease of cleaning and disinfection or sterilization ♦ Ease of handling (weight, size) and operation (ease of use, auxiliary equipment, dependence on vacuum pumps, water electricity, etc.)
Effect of the sampling device on the process or environment to be monitored	<ul style="list-style-type: none"> ♦ Disturbance of unidirectional air flow by sampling apparatus ♦ The exhaust air from the sampling apparatus should not contaminate the environment being sampled or be re-aspirated by the device ♦ Possible intrinsic addition of viable particles to the biocontamination to be measured
Efficiency and precision of the sampling method	<ul style="list-style-type: none"> ♦ Collection accuracy and efficacy
.. a d50 value smaller than 2 µm is considered appropriate ... / ...For all active microbiological air samplers, the effects of impact stress and the effect of the media dehydration during the sampling period are further considerations	<ul style="list-style-type: none"> ♦ Appropriate suction flow rate for low levels of viable airborne particles ♦ Appropriate impact/airflow velocity ♦ Sensitivity of the viable particle to the sampling procedure

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https://www.sigmaaldrich.com/technical-documents/articles/microbiology/short-webinar--the-new-en-17141-and-its-implications-for-culture.html?&utm_source=rapidmicrobiology.com&utm_medium=sponsored%20content&utm_campaign=em%20pharma-tid-2020055062

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Upcoming Webinar

ROADMAP TO COMPLIANCE -

How the new EU GMP Annex 1 (Draft 2020) is likely to impact your aseptic manufacturing

November 18, 2020
8:00 am EST, 2:00 pm CET

[Register Now](#)

Speakers:

Aram Montazami
Novatek International

Dr Anne-Grit Klees
Merck KGaA

The recent 2020 draft of EU GMP Annex 1, which regulates the manufacture, control and release of sterile pharmaceutical products in the EU, signals that there are comprehensive changes to come. This webinar will outline a holistic approach to ensure continued compliance, taking into account the expected changes, including the need for EM programs and procedures based on risk analysis.

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Comparison of guidelines and regulations of microbial levels for clean rooms											
Clean room grade ISO	Clean room grade (EU's GMP)	Clean room grade (US Fed.Std20 9E)	Active air sampling (cfu/m3)			Settle plates (diam. 90mm cfu/4 hours)			Contact plates (cfu/plate, surface 24-30 cm2)		
			EU's GMP	USP 36 <1116>	FDA 2004	EU's GMP	USP 36 <1116>	FDA 2004	EU's GMP	USP 36 <1116>	FDA 2004
Isolator / Closed RABS (ISO 5 or better)			N/A	< 0,1%	N/A	N/A	< 0,1%	N/A	N/A	< 0,1%	N/A
5	A	100 (M 3.5)	<1	< 1%	1	<1	< 1%	1	<1	< 1%	-
6	B	1000 (M 4.5)	10	< 3%	7	5	< 3%	3	5	< 3%	-
7	C	10 000 (M 5.5)	100	< 5%	10	50	< 5%	5	25	< 5%	-
8	D	100 000 (M 6.5)	200	<100	50	100	-	50	50	-	-

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
U.S. Pharmacopeia
The Standard of Quality™

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Public Health
England

Public Health England
Microbiology Services
Porton Down
Salisbury
Wiltshire
SP4 OJG

Physical and Biological Efficiency Testing of MAS-100 NT[®] Microbiological Air Sampler Using Techniques Described in ISO14698-1

Report No. 14/013 E

Report Prepared For: Merck Millipore


Issue Date: 27th January 2015

Executive Summary

The physical efficiency of the MAS-100 NT[®] air sampler for collecting small bacteria laden particles of various sizes has been compared with membrane filter samplers using the techniques described in ISO14698-1. The samplers were operated simultaneously in a controlled room where they were challenged with airborne bacteria. Uniform sized particles of different diameters containing bacterial spores were generated into the room. The results showed that the MAS-100 NT[®] sampler was approximately 60.2% efficient for collecting particles containing bacterial spores in the 0.8 micron size, in comparison with the membrane filter samplers, rising to over 90% in the range of 2 to 6 micron size. The results indicate a 50% cut-off of below 1 micron, making the sampler suitable for use within a clean room environment.

The biological efficiency of the MAS-100 NT[®] was compared to that of the Casella slit sampler, a commonly used reference sampler. The biological efficiency was measured as the comparative efficiency of collection of *Staph. epidermidis*, a common human-associated clean room contaminant, and the extremely aerostable *B. atrophaeus* spore. The biological efficiency measured was found to be 82.62% that of the Casella sampler when using the total counts, indicating the sampler is effective at sampling bacteria laden particles, without loss of viability.

Physical Efficiency of the MAS-100 NT[®] Sampler Particle Size

Report Written By

Name: Mr Simon Parks
Title: Senior Biosafety Scientist

Report Authorised By

Name: Mrs Sara Speight
Title: Senior Biosafety Scientist

Please be aware that the use of the Royal Coat of Arms is highly restricted and cannot be copied. Please do not put the PHE logo on your website or use our name to endorse your products. Any reference to PHE needs to be approved by us before it can be used.

Public Health
England

Public Health England
Microbiology Services
Porton Down
Salisbury
Wiltshire
SP4 OJG

Physical and Biological Efficiency Testing of MAS-100 ISO MH[®] (9 m Extension Tubing) Microbiological Air Sampler Using Techniques Described in ISO14698-1

Report No. 14/013 C

Report Prepared For: Merck Millipore


Issue Date: 27th January 2015

Executive Summary

The physical efficiency of the MAS-100 Iso MH[®] air sampler, with 9 metres of tubing attached, for collecting small bacteria-laden particles of various sizes has been compared with membrane filter samplers using the techniques described in ISO14698-1. The samplers were operated simultaneously in a controlled room where they were challenged with airborne bacteria. Uniform sized particles of different diameters containing bacterial spores were generated into the room. The results showed that the MAS-100 Iso MH[®] sampler was approximately 62.5 – 82.6% efficient for collecting particles containing bacterial spores in the 0.8 -1.3 micron size, in comparison with the membrane filter samplers, rising to 93.3% in the range of 5-6 micron size. The results indicate a 50% cut-off of below 1 micron, making the sampler suitable for use within a clean room environment.

The biological efficiency of the MAS-100 Iso MH[®] was compared to that of the Casella slit sampler, a commonly used reference sampler. The biological efficiency was measured as the comparative efficiency of collection of *Staph. epidermidis*, a common human-associated clean room contaminant, and the extremely aerostable *B. atrophaeus* spore. The biological efficiency measured was found to be 73.47% that of the Casella sampler when using the total counts, indicating the sampler is effective at sampling bacteria-laden particles, without loss of viability.

Physical Efficiency of the MAS-100 ISO MH[®] (with 9 meter extension) Sampler Particle Size

Report Written By

Unsigned Copy. For approval signatures please refer to original copy

Name: Mr Simon Parks
Title: Senior Biosafety Scientist

Report Authorised By

Name: Mrs Sara Speight
Title: Senior Biosafety Scientist

Please be aware that the use of the Royal Coat of Arms is highly restricted and cannot be copied. Please do not put the PHE logo on your website or use our name to endorse your products. Any reference to PHE needs to be approved by us before it can be used.

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ICR TSA+LT

- Alcohol
- Phenolics
- H₂O₂
- Peracetic acid
- Hypo chlorite
- Aldehydes
- Glucoprotamine
- QAC (approx. < 0.2%)

ICR TSA + LTHThio

- Alcohol
- Phenolics
- H₂O₂
- Peracetic acid
- Hypo chlorite
- Aldehydes
- Glucoprotamine
- Mono-Peroxyphthalate
- Dichloroisocyanurate
- QAC (approx. < 0.2%)

ICR TSA+LTHTh

- Alcohol
- Phenolics
- H₂O₂
- Peracetic acid
- Hypo chlorite
- Aldehydes
- Glucoprotamine
- Mono-Peroxyphthalate
- Dichloroisocyanurate
- QAC (approx. < 0.2%)

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
ISO 18593:2018 Microbiology of the food chain — Horizontal methods for surface sampling

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Product	Ord. No.	Shelf life in month	Slides per pack
Contact DC For disinfectant control	1.02147.0001 side A: Tryptic soy agar side B: Tryptic soy agar with neutralizer	8	10
Contact TVC For Total Colony Counts	1.02149.0001 side A: Nutrient agar with TTC side B: Nutrient agar	8	10
Contact YM(R) For the detection of Yeasts and Moulds	1.02139.0001 side A: Tryptic soy agar with TTC side B: Rose Bengal Chloramphenicol agar (R)	6	10
Contact E For the detection of Enterobacteriaceae	1.02137.0001 side A: Plate count agar side B: VRBD agar	6	10
Contact C For the detection of Total coliforms/E.coli	1.02136.0001 side A: Plate count agar side B: Chromocult® Coliform agar	8	10

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