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

ISO 8199:2018
Water quality -- General requirements and guidance for microbiological examinations by culture

8.2.5 Membrane filtration method

8.2.5.1 Test portion

The maximum volume of the test portion depends on the filterability of the water sample and on the membranes used. This technique is suitable for waters which contain little particulate matter in suspension, for example drinking water. With membranes of mean pore diameter 0,45 μm , it may be possible

3. ENUMERATION METHODS

Use the membrane filtration method or the plate-count methods, as prescribed. The most-probable-number (MPN) method is generally the least accurate method for microbial counts, however, for certain product groups with a very low bioburden, it may be the most appropriate method.

c. Membrane filter method: The membrane filter method (9215D) can be used to analyze large volumes of low-turbidity water and is the method of choice for waters with low numbers of heterotrophic organisms (<1 to 10 CFU/mL). It produces no heat shock. However, there is the expense of the membrane filter, the smaller area for colony development, possible damage to cells by excessive filtration pressures, possible variations in membrane-filter quality (see Section 9020B.5i), and the need to detect colonies via reflected light against a white background if colored filters or contrast stains are not used.

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S-PAK® MEMBRANE Certificate of Quality

We certify that the product described within meets the following criteria.

Catalogue Number : HAWG047S6

Lot Number : F3PA84602

Expiry date : NOV-2015

Product Description

The membrane used in this product (0.45 µm mixed esters of cellulose) has been sampled and tested and conforms to the specifications.

The membrane filter has been manufactured for use in the bacteriological analysis of potable waters in accordance with procedures referenced in Standard method (recent edition), and laboratory standards for equipment and materials set forth by the U.S. Environmental Protection Agency. Each lot is tested for recovery in accordance with ISO 7704 procedure and meets or exceeds the highest suggested recoveries.

Good Manufacturing Practices

This product was manufactured in a Millipore SAS facility which adheres to Good Manufacturing Practices.

ISO® 9001 Quality Standard

This product was manufactured in a Millipore SAS facility whose Quality Management System is approved by an accredited registering body to ISO 9001 Quality System Standard.

HAWG0--S6C0Q Rev.E 01/13
S-Pak® is a Registered Trademark of Merck KGaA
Milli Q® is a Registered Trademark of Merck KGaA
ISO® is a Registered Trademark of the International Organization for Standardization
ATCC® is a Registered Trademark of the American Type Culture Collection
The M Mark is a trademark of Merck KGaA, Darmstadt, Germany.
©2012 Merck KGaA, Darmstadt, Germany. All rights reserved.

ISO® 14001 Environmental Management

This product was manufactured in a Millipore SAS facility whose Environmental Management program is approved by an accredited registering body to the appropriate ISO 14001 System Standard.

Quality Assurance Lot Release Criteria

This manufacturing lot was sampled, test and released by Quality Assurance for compliance with the following characteristics :

Statistical Controls

- Physicals Tests
- Use tests of the filtration device
- Conformity of labeling and packaging.

Biological tests

Sterility

This product has been sterilized by Ethylene Oxide (EO) in a validated sterilization cycle. Biological indicators incorporated in the lot have shown to be sterile.

Fertility / Growth promotion

Criteria	Specifications	Results
<i>Escherichia coli</i> ATCC® 8739	Recovery ≥ 90% versus spread plate	Conform
<i>Enterobacter aerogenes</i> ATCC® 49701	Recovery ≥ 90% versus spread plate	Conform

Lot Analysis

This product was designed and manufactured to meet the following specifications. Performance is confirmed by testing on each lot.

Membrane Filter Analysis

- Tested on a lot audit basis :

Criteria	Specifications	Results
Pore size	0.45 µm (Determined by retention of <i>Serratia marcescens</i>)	Conform
Water Extractables	≤ 1.0%	Conform
Membrane thickness	115 – 180 µm	Conform
Water flow time of the membrane	25-50 sec/500ml @47mm at 27.5 ± 0.5 inches Hg	Conform
Wettability of the membrane	5 s maximum with Milli Q® water	Conform
Extractable Inorganic Na	≤ 0.01 mg/filter	Conform

- The manufacturing process has been validated to conform with the following specifications :

Criteria	Specifications	Results
Extractables Inorganics :		
Ca	≤ 0.01 mg/filter	Conform
Zn	≤ 0.001 mg/filter	Conform
Ag	≤ 0.0001 mg/filter	Conform
Pb	≤ 0.0001 mg/filter	Conform
Hg	≤ 0.0001 mg/filter	Conform
Cu	≤ 0.0001 mg/filter	Conform
Cr	≤ 0.0001 mg/filter	Conform
Organic (TOC)	≤ 1.0 mg/filter	Conform
P (reported as Phosphate)	≤ 0.02 mg/filter	Conform
Ammonia (NH ₃)	≤ 0.01 mg/filter	Conform
Neutrality	pH of extract ≤ 8.3	Conform
Total acidity	≤ 0.001 mEq/filter	Conform

According to the above results, the product complies with Millipore SAS's acceptance criteria and is released.

Christel Noehringer
BioMonitoring Quality Manager

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0,45 µm
MCE
100 & 250ml

Water testing, Bioburden and Pharma water, F&B application with liquid media (presence of pad to convert the unit into a Petri dish)

0,22 µm
MCE
100ml

Water systems or water loops validation, product bioburden testing, pharma water microbial testing

0,45 µm
MCE
100 & 250 ml

Water testing, Bioburden and Pharma water, Hard to filter turbid matrix, price sensitive pharma customers (multipack)

0,45 µm
PVDF
100 ml

Low binding membrane for Pharma bioburden, antibiotics and solvents testing, or samples not chemically compatible with MCE

0,8 µm
MCE
100 & 250 ml

Food & Beverages yeast and moulds testing : wine, beer, cider, soft drinks

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MILLIFLEX® PRE-FILLED MEDIA CASSETTES

Catalogue Number : MXSMCTS48
 Medium Abbreviation : TRYPTIC SOY
 Lot Number :
 Expiration Date :

Product Description

Each cassette is a plastic 3 piece assembly containing approximately 10ml of Tryptic Soy Agar.

Composition of the Medium

Formula in grams per liter of purified water :

Peptone from Casein	15.0 g
Peptone from Soymeal	5.0 g
Sodium Chloride	5.0 g
Agar	15.0 g

This medium is recommended for a variety of microorganisms.
 Storage conditions: 2 - 8°C

Good Manufacturing Practices

This product was manufactured in a Millipore facility which adheres to Good Manufacturing Practices.

ISO® 9001 Quality Standard

This product was manufactured in a Millipore facility whose Quality Management System is approved by an accredited registering body to ISO 9001 Quality System Standard.

Pharmacopoeia

This product was manufactured, tested and released according to European and US Pharmacopoeia recommendations.

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This manufacturing lot was sampled, tested and released according to the following specifications :

100 % Inspection

Visual control of appearance and medium level of each cassette.

Statistical Controls

- * Packaging integrity
- * Conformity of labelling and packaging
- * Use test of the cassette.

Biological Tests

Sterility

Representative samples were subjected to sterility test by direct incubation of the pre-filled cassettes.

Fertility / Growth Promotion

Fertility tests were conducted by direct inoculation of the media. All the samples provided good growth and typical colony morphology when challenged with approximately 10 - 100 CFU per cassette.

Batch record

The lot of cassettes has a complete documentation which includes product description, test protocol, test specification and results.

Lot Analysis

This product was designed and manufactured to meet the following specifications.

Criteria	Specifications	Results
Appearance of the media	light amber, slightly opalescent	Conforms
pH	7.3 +/- 0.2	Conforms
Weight of Agar	Min : 10.0 g Max : 12.0 g	Conforms
Sterility Assurance Level	≤ 1 % contamination widest after 7 days at 20-25°C and 30-35°C	Conforms
Growth promotion test		
Staphylococcus aureus ATCC® 6538	Good growth after 24 h at 30-35°C	Conforms
Bacillus subtilis ATCC® 6633	Good growth after 24h at 30-35°C	Conforms
Pseudomonas aeruginosa ATCC® 9027	Good growth after 24h at 30-35°C	Conforms
Escherichia coli ATCC® 8739	Good growth after 24h at 30-35°C	Conforms
Candida albicans ATCC® 10231	Good growth after 48h at 20-25°C	Conforms
Aspergillus niger ATCC® 16404	Good growth after 5 days at 20-25°C	Conforms

According to the above results, the product complies with Millipore's acceptance criteria and is released.

Brigitte Weber
 Quality Assurance Manager

PF08952, Rev.
 Millipore and Milliflex are registered trademarks of Millipore Corporation.
 ISO is a registered trademark of the International Organization for Standardization or an affiliated company.
 ATCC is a registered trademark of American Type Culture Collection.

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Will end
New

Hardware, accessories, spare parts and services	Milliflex	Milliflex Oasis
Filtration pump, one system with ONE filtration heads	MXPPLUS01	
Filtration pump, one system with TWO filtration heads	MXPPLUS02	MMSYSTEMM1
Filtration pump, one system with THREE filtration heads	MXPPLUS03	
Filtration pump, with TWO filtration heads for EZ-Fit Filtration Units		MMSYSTEMZ1
Filtration head adapter for Milliflex PLUS pump for Mfx Oasis funnels		MXHEADMM1

Filtration units	Milliflex	Milliflex Oasis
100 mL funnel, 0.45 µm white gridded, mixed cellulose esters (MCE) membrane	MXHAWG124	MMHAWG124
250 mL funnel, 0.45 µm white gridded, mixed cellulose esters (MCE) membrane	MXHAWG224	MMHAWG224
100 mL funnel, 0.45 µm black gridded, mixed cellulose esters (MCE) membrane	MXHABG124	MMHABG124
100 mL funnel, 0.22 µm white gridded, mixed cellulose esters (MCE) membrane	MXGSWG124	MMGSWG124
100 mL funnel, 0.45 µm white plain PVDF Durapore® membrane	MXHVWP124	MMHVWP124
250 mL funnel, 0.45 µm white plain, PVDF Durapore® membrane	MXHVWP224	MMHVWP224

Culture media	Milliflex	Milliflex Oasis
R2A low nutrient media plates	MXSMCRA48	MMSMCRA48
Tryptic soy media plates for detection of aerobic and anaerobic flora	MXSMCTS48	MMSMCTS48
Sabouraud dextrose media plates for detection of yeast & molds	MXSMCSD48	MMSMCSD48
Milliflex® cassette sabouraud dextrose agar (w/Chloramphenicol)	MXSMCSP48	
Milliflex® cassette m-Endo agar for the detection of coliforms	MXSMEND48	
Milliflex® cassette PCA agar for the detection of bacteria in water	MXSMPCA48	
Milliflex® cassette cetrimide agar with naladixic acid	MXSMCET48	
Milliflex® cassette TSA and lecithin agar	MXSMPLP48	
Milliflex® cassette HPC agar media, for detection of plate count bacteria	MXSMHPC48	
Milliflex® cassette KF strep agar for the detection of fecal streptococci	MXSMKFS48	
Milliflex® cassette PIA agar for the detection and recovery of P aeruginosa	MXSMPIA48	
Mac Conkey Milliflex® agar for detection of lactose fermenting bacteria	MXSMCMC24	
Milliflex® empty cassette for agar media	MXSMC0120	
Milliflex® empty cassette for liquid media	MXLMC0120	

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Hardware, accessories, spare parts and services	Milliflex	Milliflex Oasis
Filtration pump, one system with ONE filtration heads	MXPPLUS01	
Filtration pump, one system with TWO filtration heads	MXPPLUS02	MMSYSTMM1
Filtration pump, one system with THREE filtration heads	MXPPLUS03	
Filtration pump, with TWO filtration heads for EZ-Fit Filtration Units		MMSYSTEME1
Filtration head adapter for Milliflex PLUS pump for Mfx Oasis funnels		MXHEADMM1

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Filtration units	Item #	END DATE
Milliflex® 100 mL funnel, 0.45 µm white gridded, MCE membrane	MXHAWG124	Q1 2022
Milliflex® 250 mL funnel, 0.45 µm white gridded, MCE membrane	MXHAWG224	Q2 2022
Milliflex® 100 mL funnel, 0.45 µm black gridded, MCE membrane	MXHABG124	Q1 2022
Milliflex® 100 mL funnel, 0.22 µm white gridded, MCE membrane	MXGSWG124	Q1 2022
Milliflex® 100 mL funnel, 0.45 µm white plain PVDF Durapore® membrane	MXHVWP124	Q2 2022
Milliflex® 250 mL funnel, 0.45 µm white plain, PVDF Durapore® membrane	MXHVWP224	Q2 2022
Culture media		
Milliflex® cassette R2A low nutrient media	MXSMCRA48	Q1 2022
Milliflex® cassette tryptic soy media plates for detection of aerobic and anaerobic flora	MXSMCTS48	Q1 2022
Milliflex® cassette sabouraud dextrose media plates for detection of yeast & molds	MXSMCSD48	Q1 2022
Hardware, accessories, spare parts and services		
Milliflex PLUS pump single head	MXPPLUS01	Q1 2021
Milliflex PLUS pump double head	MXPPLUS02	Q2 2020
Milliflex PLUS pump triple head	MXPPLUS03	Q2 2020
Spare parts for Milliflex PLUS		Q1 2024

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Filtration units	Item #	AVAILABILITY
Milliflex Oasis® 100 mL funnel, 0.45 µm white gridded, mixed cellulose esters (MCE) membrane	MMHAWG124	Launch date
Milliflex Oasis® 250 mL funnel, 0.45 µm white gridded, mixed cellulose esters (MCE) membrane	MMHAWG224	June 2020
Milliflex Oasis® 100 mL funnel, 0.45 µm black gridded, mixed cellulose esters (MCE) membrane	MMHABG124	March 2020
Milliflex Oasis® 100 mL funnel, 0.22 µm white gridded, mixed cellulose esters (MCE) membrane	MMGSWG124	March 2020
Milliflex Oasis® 100 mL funnel, 0.45 µm white plain PVDF Durapore® membrane	MMHVWP124	June 2020
Milliflex Oasis® 250 mL funnel, 0.45 µm white plain, PVDF Durapore® membrane	MMHVWP224	June 2020
Milliflex Oasis® rapid 100 mL funnel, 0.45 µm white plain, PVDF Durapore® membrane	MMHVAFX24	June 2020
Culture media		
Milliflex Oasis® R2A low nutrient media plates	MMSMCRA48	Launch date
Milliflex Oasis® tryptic soy media plates for detection of aerobic and anaerobic flora	MMSMCTS48	Launch date
Milliflex Oasis® sabouraud dextrose media plates for detection of yeast & molds	MMSMCSD48	Launch date
Hardware, accessories, spare parts and services		
Milliflex Oasis® filtration pump, one system with two filtration heads and no power supply	MMSYSTEMM1	Available
Filtration head adapter for Milliflex PLUS pump for Mfx Oasis funnels	MXHEADMM1	March 2020

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Millipore.
Preparation, Separation,
Filtration & Testing Products

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“For microbial recovery and identification, microbiological testing laboratories sometimes use alternative test methods to those described in the general chapters for a variety of reasons, such as economics, throughput, and convenience. **Validation of these methods is required**” (EP, USP)

Contents

Method Summary	
Introduction	7
Milliflex Quantum Method Description	8
References	9
Test Summaries	10
Materials	12
Methods	13
Glossary	17
Millipore Qualification Test Methods & Results	
Robustness	21
Ruggedness	24
Accuracy	28
Precision	30
Linearity	31
Range	32
Limit of Quantification	33
Limit of Detection	34
Specificity	36
False Positives Study	37
Conclusion	38

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MERCK

Milliflex Quantum Rapid Detection System Validation Protocol
MXQUA4VP1 – Rev A
Installation Qualification

Milliflex Quantum Rapid Detection System Validation Protocol
MXQUA4VP1 – Rev A
Installation Qualification

INSTALLATION QUALIFICATION PROTOCOL FOR THE MILLIFLEX® QUANTUM RAPID DETECTION SYSTEM

Table of Contents

INSTALLATION QUALIFICATION PROTOCOL FOR THE MILLIFLEX® QUANTUM RAPID DETECTION SYSTEM	36
1. OBJECTIVE	39
2. SCOPE	39
3. RESPONSIBILITIES	39
4. SAFETY	39
4.1 Protective Equipment	39
4.2 Incident Documentation	39
5. PREREQUISITES	40
5.1 Reference Documents	40
5.2 Reference Instruments	40
5.3 Purchasing Information	40
5.4 Utilities and Operating Environment Requirements	40
5.5 Work Area, Equipment and Personnel Preparation	41
5.6 Test Material Requirement	42
6. TEST PROCEDURES	43
6.1 Verification and Identification of the Milliflex Quantum Rapid Detection System	43
6.2 Verification of the Milliflex Quantum Membrane Transfer Tool Assembly	45
6.3 Verification of the Milliflex Quantum Reader Assembly and Start-Up	47
6.4 Verification of the Milliflex Quantum Spot Counter Software Installation, Camera Assembly and Start-Up	51
6.5 Verification of the Milliflex Quantum Spot Counter Software Version Installed	54
6.6 Installation Qualification Summary	57
Appendix 1 Signature Identification Sheet	58
Appendix 2 Reference Documents	59
Appendix 3 Reference Instruments	60
Appendix 4 Purchasing Information	61
Appendix 5 Utilities and Operating Environment Requirements	62
Appendix 6 Work Area, Equipment and Personnel Preparation	63
Appendix 7 Verification and Identification of the Milliflex Quantum Rapid Detection System	64
Appendix 8 Verification of the Milliflex Quantum Membrane Transfer Tool Assembly	67
Appendix 9 Verification of the Milliflex Quantum Reader Assembly and Start-Up	68
Appendix 10 Verification of the Milliflex Quantum Spot Counter Software Installation, Camera Assembly and Start-Up	70
Appendix 11 Verification of the Milliflex Quantum Spot Counter Software Version Installed	72
Appendix 12 Installation Qualification Summary	73
Appendix 13 Comment Sheet	74

Original
is
blue

Copyright 2010 by Millipore Corporation.
Copying of this page for any purpose is a violation of copyright laws.
← Original (by color: Blue = Original page/ Black = Copy

Page 36
of 219

Original
is
blue

Copyright 2010 by Millipore Corporation.
Copying of this page for any purpose is a violation of copyright laws.
← Original (by color: Blue = Original page/ Black = Copy

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MERCK

PDA Journal
of Pharmaceutical Science and Technology

Validation of Milliflex® Quantum for Bioburden Testing of Pharmaceutical Products

Oliver Gordon, Marcel Goverde, Alexandra Staerk, et al.

PDA J Pharm Sci and Tech 2017, 71 206-224
Access the most recent version at doi:[10.5731/pdajpst.2016.007450](https://doi.org/10.5731/pdajpst.2016.007450)

ABSTRACT: This article reports the validation strategy used to demonstrate that the Milliflex® Quantum yielded non-inferior results to the traditional bioburden method. It was validated according to USP <1223>, European Pharmacopoeia 5.1.6, and Parenteral Drug Association Technical Report No. 33 and comprised the validation parameters *robustness, ruggedness, repeatability, specificity, limit of detection and quantification, accuracy, precision, linearity, range, and equivalence in routine operation*. For the validation, a combination of pharmacopeial ATCC strains as well as a broad selection of in-house isolates were used. In-house isolates were used in stressed state. Results were statistically evaluated regarding the pharmacopeial acceptance criterion of $\geq 70\%$ recovery compared to the traditional method. Post-hoc test power calculations verified the appropriateness of the used sample size to detect such a difference. Furthermore, equivalence tests verified non-inferiority of the rapid method as compared to the traditional method. In conclusion, the rapid bioburden on basis of the Milliflex® Quantum was successfully validated as alternative method to the traditional bioburden test.

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