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
MC-Media Pads

E.Coli & Coliforms

In order to showcase the capabilities of our product, we have previously submitted a new accreditation under AOAC and ISO which covers a wider range of Food matrices. Since the completion of the task, we have been requested by the regulatory body to alter the article information of our product. the primary goal was to avoid confusion with our accreditation, new and previous.

Also, with the change of article code, we took this opportunity to revalidate the shelf life. Thus, the shelf life is changed from 36 to 18 months. This will be returned to 36 months, again, once the revalidation is completed.

Please be informed that the current article code 1323570001 MC-Media Pad EC will be available until stocks exist. Once depleted, we will transfer all orders to the new article code 1323570002.

Product Description	Article Code	Accreditation	Shelf life
MC-Medi Pad EC	1.32357.0001	AOAC® Performance Tested SM certificate number 070901	36 months
	1.32357.0002	AOAC® Performance Tested SM certificate number 011902 Microval Summary report 2017LR71	18 months

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Millipore®

Preparation, Separation,
Filtration & Testing Products

MC-Media Pad

How to handle difficult to test matrices

The MC-Media Pad is designed as a convenient method for rapid routine testing of microbial contamination in food and beverage samples throughout production, from raw materials, to finished products. The MC-Media Pad solution is composed of a series of ready-to-use pads for total count and specific detection and enumeration of indicator organisms.

Each MC-Media Pad consists of a fabric pad coated with a dedicated culture media formulation, placed on an adhesive support with a transparent gas permeable cover.

MC-Media Pads have been specifically developed for the testing of 1 mL food samples.

An extensive list of food and beverage matrices described in the AOAC-PTM and MicroVal validation reports¹ and detailed in the document entitled MC-Media Pad: A Convenient and Approved Method for Counting Contaminants in Food and Beverage products², have been tested with the MC-Media Pad product range to demonstrate equivalency with control methods. The

Application Note

Conclusion:

Results generated in this study demonstrate that MC-Media Pads can be used to easily detect and enumerate indicator organism contamination in challenging food and beverage matrices with properties including high sugar or fat content, high viscosity, high acidity or enzymatic or colored backgrounds.

In addition to the suitability for difficult to test matrices, MC-Media Pads provide some benefits in terms of color-coding, handling, time to results for yeasts and molds in 48 hours, and improved read-out using chromogenic substances / specific dyes allowing a fast and convenient control of microbial contamination.

Food products	Characteristic
Mozzarella	Cheese products with potential enzymatic activity
Brie	
Powdered Pecorino	
Flavored yogurt	Dairy product
Black pepper	Dehydrated product
Dry soup	
Chocolate pudding	Cooked product
Cola	Acidic and high sugar content products with a black colored background
Inverted sugar syrup	Viscous product with high sugar content
Mayonnaise	Fatty and viscous product

Sample Preparation

Portions of 10 g or 10 mL of each food sample were added to a stomacher bag or blender cup. 90 mL of

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Millipore®

Preparation, Separation,
Filtration & Monitoring Products

Suitability of MC-Media Pad® for Hygiene Monitoring

(passive air, personnel and surface monitoring)

This study shows that the MC Media Pad® is suitable for a range of hygiene monitoring applications.





The study focused on usage of MC-Media Pad® in passive air, personnel and surface monitoring in comparison to standard test methods.

The MC-Media Pad® is designed as a convenient method for rapid routine testing of microbial contamination in food and beverage samples throughout

production, from raw materials to finished products. The MC-Media Pad® solution is composed of a series of ready-to-use pads for total count and specific detection and enumeration of indicator organisms.

Each MC-Media Pad® consists of a fabric pad coated with a dedicated culture media formulation, placed on an adhesive support with a transparent gas-permeable cover (figure 1). The MC-Media Pad® has been specifically developed for the testing of 1 mL samples.

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