

Revolution in Cell Culture from **biowest**



biowest) Instant FBS

Instant FBS is a powdered, soluble Foetal Bovine Serum offering many unique benefits :

Easy to use :

- Dilution in cell culture water or directly in cell culture Media
- 3 g of instant FBS is equivalent to 50 ml of liquid FBS
- 3 g prepacked in tubes for 500 ml Media (10% FBS)

Quality :

- Same quality as liquid FBS
- No more aliquoting and refreezing of FBS
- Minimized risk of contamination

Storage :

- Instant FBS is stored at 2-8°C for a period of 3 years.

Environmental :

- No more shipment on dry ice, reducing cost and CO₂ pollution
- No storage in energy consuming freezers in the lab
- Less weight and volume transported ; no freezer truck, for reduced transport cost and less environmental contamination



Revolution in Cell Culture from **biowest**



Request for free Sample at Fax Nr. 0033 2 41 46 38 90 :

Contact Name

Company :

Address

Tel. / Fax :

E-Mail :

We want to test Instant Sterile FBS :

Instant Sterile FBS S381G-N0.05L (for 50 ml)

Quantity of samples requested :

Expected consumption :

You need further informations concerning our product line ? Don't hesitate to contact us !



Biowest S.A.S. Rue de la Caille | 49340 Nuillé | France
Tél. : [+33] 241 464 242 | Fax : [+33] 241 463 890 | Email : biowest@biowest.net