

PACK
5
SACHETS

UNIQUE
PRESENTATION



**EXACTLY
WHAT YOU NEED.**

**SACHETS OF PREWEIGHED DEHYDRATED
CULTURE MEDIA SUFFICIENT TO PREPARE 500ml.**

Scharlau
The wise choice

QUICK

Sachets contain the exact amount of media to be dissolved in 500mL of water, eliminating weighing errors and saving time.

SAFE


The handling of large quantities of culture media is no longer necessary and the risk of possible inhalation and exposure to dust is eliminated.

DURABLE

Being impermeable, sachets are resistant to moisture, and product integrity is maintained , providing a shelf life of 4 to 5 years.

ECONOMIC


Presented in boxes of 5 sachets, this packaging makes the product ideal when large quantities of media are not required. In this way, you only pay for the exact amount that you need .



See video tutorial on our website: www.scharlab.com

Avoid weighing errors.

The exact quantity that you will use.



No waste created.

Ideal for unexpected peaks in workflow.

Ordering information:

Description	Art. No	Description	Art. No
Acetate Differential Agar	01-665BA05	Microinstant® Chromogenic	
Antibiotic Medium A pH		Listeria Agar acc. to Ottaviani & Agosti	01-719BA05
6,6 EUR.PHARM.	01-009BA05	Microinstant® Sakazakii	
Antibiotic Medium A pH 7,9 eur. pharm.	01-017BA05	Chromogenic Agar	01-744BA05
Asparagine Broth	02-271BA05	Microinstant® Tryptone	
Bacillus Cereus Agar	01-262BA05	Bile Glucuronid Agar (TBX Agar)	01-619BA05
Bacillus Cereus Selective Agar	01-487BA05	Mineral Modified Glutamate Agar	01-659BA05
Baird Parker Agar	01-030BA05	Mineral Modified Glutamate Medium	02-656BA05
BAT (Bacillus acidoterrestris) Agar	01-675BA05	Motility Nitrate Medium	03-612BA05
BAT (Bacillus acidoterrestris) Broth	02-675BA05	MRS Agar	01-135BA05
BCYE Legionella Agar	01-687BA05	MRS Broth	02-135BA05
Beerens Cosmetic Diluent	02-257BA05	MRVP Broth	
BGA (Brilliant Green Agar)	01-203BA05	(Methyl Red Voges Proskauer Broth)	02-207BA05
BGA Modified		Mueller-Hinton Agar	01-136BA05
(Brilliant Green Modified Agar)	01-309BA05	Mueller-Hinton Broth	02-136BA05
BGBL Broth		Müller-Kauffmann Tetrathionate Broth	02-335BA05
(Brilliant Green Bile 2% Broth)	02-041BA05	Neutralizing Fluid	02-512BA05
BHI Agar (Brain Heart infusion Agar)	01-599BA05	Nickerson Agar (BIGGY Agar)	01-137BA05
BHI Broth (Brain Heart infusion Broth)	02-599BA05	Nutrient Agar	01-140BA05
Bile Esculin Azide Agar	01-592BA05	Nutrient Agar (APHA)	01-144BA05
Bile Esculin Modified Agar	01-265BA05	Nutrient Broth	02-140BA05
Blood Agar	01-352BA05	Nutrient Broth (APHA)	02-144BA05
Blood Agar (Columbia)	01-034BA05	Oxford Agar	01-471BA05
Blood Agar No. 2	01-505BA05	Palcam Agar	01-470BA05
Bolton Enrichment Broth	02-688BA05	Phenylalanine Agar	01-083BA05
Buffered Peptone Water	02-277BA05	Plate Count Agar (PCA)	01-161BA05
Buffered Peptone Water Eur. Pharm.	02-494BA05	Plate Count Modified Agar	01-329BA05
Burkholderia cepacia		Plate Count Skim Milk Agar	01-412BA05
Selective Agar BCSA	01-736BA05	Potato Dextrose Agar (PDA)	01-483BA05
Cary-Blair Transport Medium	03-643BA05	Potato Dextrose Broth (PDB)	02-483BA05
Casein Lecithin Polysorbate Broth	02-539BA05	R2A Agar	01-540BA05
Cetrimide Agar	01-160BA05	Rappaport Vassiliadis Broth	02-379BA05
CGA (Chloramphenicol Glucose Agar)	01-366BA05	Rappaport Vassiliadis	
Chapman TTC Agar	01-053BA05	Modified Semi-Solid Medium	03-376BA05
CLED Agar	01-047BA05	Rappaport Vassiliadis	
CN (Cetrimide-Nalidixic)		Salmonella Enrichment Broth	02-668BA05
Selective Agar (01-609)	01-609BA05	Reinforced Clostridial Medium (RCM)	03-289BA05
Columbia Agar (Eur. Pharm.)	01-680BA05	Rose Bengal Agar	01-301BA05
D/E Neutralizing Broth	02-610BA05	Sabouraud Broth	02-165BA05
Deoxycholate Lactose Agar	01-057BA05	Sabouraud Chloramphenicol Agar	01-166BA05
Dermatophyte Selective Agar	01-442BA05	Sabouraud Dextrose Agar	01-165BA05
Dichloran Rose Bengal Chloramphenicol Agar (DRBC Agar)	01-657BA05	Sabouraud Oxytetracycline	
DNase Agar	01-346BA05	Agar (OGYEA)	01-275BA05
E. E. Broth		Salmonella Shigella Agar (SS Agar)	01-555BA05
(Enrichment Enterobacteriaceae Broth)	02-064BA05	Selenite Broth	02-598BA05
EC Broth (E. coli Broth)	02-060BA05	Selenite Cystine Broth	02-602BA05
EMB Agar (Eosin Methylene Blue Agar)	01-068BA05	Shigella Broth	02-662BA05
Endo Agar	01-589BA05	SIM Medium	03-176BA05
EUGON LT 100 Agar	01-654BA05	Simmons Citrate Agar	01-177BA05
EUGON LT 100 Broth	02-654BA05	Slanetz Bartley Agar	01-579BA05
Fraser Listeria Enrichment Broth	02-496BA05	SPS Agar	01-050BA05
Giolitti-Cantoni Broth	02-230BA05	Stuart Ringertz Transport Medium	03-454BA05
Glucose Bromocresol Purple Agar	01-502BA05	TBA (Tryptose Bile Agar)	01-526BA05
Gram Negative Broth (GN Broth)	02-093BA05	TCBS Agar	01-567BA05
Hektoen enteric Agar	01-216BA05	Tetrathionate Bile Brilliant green Broth (EUR. Pharm.)	02-629BA05
Iron Sulfite Modified Agar	01-634BA05	Thioglycollate Broth	02-186BA05
K Agar	01-674BA05	Thioglycollate Fluid Medium	03-187BA05
King A Agar (P Agar)	01-001BA05	Tryptone Water	03-156BA05
King B Agar (F Agar)	01-029BA05	Tryptone Yeast Extract Agar	01-590BA05
Kligler Iron Agar (KIA)	01-103BA05	Tryptose LaurylSulfate Broth	02-108BA05
Lactose Broth	02-105BA05	TSA (Tryptic Soy Agar)	01-200BA05
LB Agar according to Miller	01-385BA05	TSA Lecithin & Polysorbate	01-613BA05
LB Broth according to Miller	02-385BA05	TSB (Tryptic Soy Broth)	02-200BA05
Lethen Broth	02-236BA05	TSB w/o Dextrose	02-227BA05
Lethen Modified Agar	01-237BA05	TSC (Tryptose Sulfite Cycloserine Agar)	01-278BA05
Lethen Modified Broth	02-237BA05	TSI (Triple Sugar Iron Agar)	01-192BA05
Lovett Listeria Enrichment Broth	02-498BA05	TSN (Tryptone Sulfite Neomycin Agar)	01-195BA05
Lysine Iron Agar (LIA)	01-094BA05	Urea Agar	01-261BA05
m-CP Agar (Clostridium perfringens)	01-513BA05	Urea Broth	02-202BA05
m-Green Broth	02-633BA05	UVM Listeria Enrichment Broth	02-472BA05
MacConkey Agar	01-118BA05	Vogel Johnson Agar	01-206BA05
MacConkey Broth	02-118BA05	VRB Agar (Violet Red Bile Agar)	01-164BA05
MacConkey Broth (EUR. PHARM.)	02-611BA05	VRBD (Violet Red Bile Dextrose Agar)	01-295BA05
Malt Extract Agar No. 1	01-111BA05	XLD Agar	
Malt Extract Agar No. 2	01-573BA05	(Xylose Lysine Deoxycholate Agar)	01-211BA05
Mannitol Salt Agar	01-116BA05	XLD Modified Agar (Xylose Lysine	
Maximum Recovery Diluent (MRD)	02-510BA05	Deoxycholate Modified Agar)	01-552BA05
Microinstant® Chromogenic		Yersinia CIN Agar	01-444BA05
Coliform Agar Base	01-695BA05	YSG Agar (Yeast Starch Glucose Agar)	01-673BA05
Microinstant® Chromogenic			
Colinstant Agar	01-618BA05		

Scharlau is the division of **Scharlab, S.L.** and is the brand which encompasses chemicals, microbiology, glassware and chromatography. We have manufactured laboratory reagents since 1949 and microbiology culture media since 1979, being represented in more than 100 countries. We can offer you a solution for any laboratory chemical, microbiology or consumable challenges that you may have.



Scharlab S.L.

Gato Pérez, 33. Pol. Ind. Mas d'en Cisa.
08181 Sentmenat, Barcelona
Tel.: +34 93 715 19 40 - Fax: +34 93 715 27 65
E-mail: export@scharlab.com

Scharlab Italia S.r.l.

Via Alcide De Gasperi 56.
20070 Riozzo Di Cerro al Lambro (Mi), Italy
Tel.: +39 02 9823 0679 / +39 02 9823 6266
Fax: +39 02 9823 0211 / +39 02 9811 9288
E-mail: customerservice@scharlab.it

Scharlab Magyarország Kft.

4034 Debrecen, Vágóhíd. u. 2., Hungary
Tel: 0036-787-634 - Fax: 0036-781-081
E-mail: info@scharlab.hu

Scharlab Philippines, Inc.

18G Miller Compound, Barangay Bungad,
San Francisco Del Monte, Quezon City
1105 Philippines
Tel.: +63 2 3514972 - Fax: +63 2 3514972
E-mail: infophilippines@scharlab.ph

Scharlau Trading (Hangzhou) Company Ltd.

Room 508, Modern Star, 337# Shaoxing Road,
Hangzhou, Zhejiang, China
Tel: 0086-571-88536628
Fax: 0086-571-88536698
E-mail: china01@scharlab.com

www.scharlab.com