

FBS, Qualified

Lot Number: 2440087
Item Number: 10270
Expiration Date: MAY2027
Storage Temp: <= -10C
Country of Origin: Brazil

TEST	TEST ID	SPECIFICATION	RESULT	UNITS
Albumin	ELECTR0029	1.70 - 3.50	2.38	g/100ml
Alpha Globulin	ELECTR0030	0.70 - 2.00	1.41	g/100ml
Appearance (Liquid)	SOLAPP0066	Satisfactory	Satisfactory	
Bacterial	STERI0024	Negative	Negative	
Beta Globulin	ELECTR0032	0.30 - 0.90	0.33	g/100ml
BVD	VIRUS0159	Tested	Tested	
Diploid Growth	HDFG0001	0.70 - 1.40	0.90	
Documentation	DOC0014	Satisfactory	Satisfactory	
Electro Profile	ELECTR0035	Typical	Typical	
Endotoxin	ENDO0068	0.000 - 50	0.202	EU/mL
Fungal	STERI0027	Negative	Negative	
Gamma Globulin	ELECTR0034	10.00 - 200.00	38.55	ug/ml
Haemoglobin	HAEMO0023	1.00 - 30.00	12.29	mg/100ml
Health Certificate Number	MISC0493	Check and Record	00671/CC-MET/22	
IBR	VIRUS0161	Negative	Negative	
Mycoplasma	MYCO0026	Negative	Negative	
Osmolality	OSMO0025	280 - 340	312	mOsm/kg H2O
pH	PH0037	6.90 - 7.80	7.43	pH Units
PI3	VIRUS0160	Negative	Negative	
Relative Cloning Efficiency	CLON0003	0.70 - 1.40	1.00	
Relative Myeloma/Hybridoma	CCMISC0021	0.70 - 1.40	0.82	
Relative Plating Efficiency	PLATE0006	0.70 - 1.40	0.95	
Supplier Lot Number	MISC0492	Check and Record	BIO/22/002	
Total Protein	PROTEI0035	3.50 - 5.50	4.12	g/100ml
Toxicity	CYTTOX0008	Satisfactory	Satisfactory	

Quality Systems Department

Date: 13-Jun-2022

References

- DOC0014: BMR review carried out - Test Reference